




Distal Radius Plate Instrument Set

Sr	Product No.	Name	Specification	QTY	Unit	Outview
1	09S7001	Guide Pin	$\phi 1.5 \times 150$	2	pcs	
2	09S7002	Thread Pin	$\phi 1.5 \times 150$	2	pcs	
3	09S7003	Pin Sleeve	$\phi 1.5$	2	pcs	

4 09S7004 Drill Bit $\phi 1.8$ 3 pcs



5 09S7005 Drill Sleeve $\phi 1.8$ 2 pcs



6 09S7006 Drill Bit $\phi 2.0$ 2 pcs



7 09S7007 Drill Sleeve $\phi 2.0 - \phi 2.7$ 1 pcs



8	09S7008	Dril Bit	$\phi 2.7$	1	pcs
9	09S7009	Tap	HA2.7	1	pcs
10	09S7010	Stardrive		1	pcs
11	09S7011	Instrument Box		1	set
		In total		18	pcs



Distal Radius Plate Instrument Set

