

## LCP Small Fragment Instrument Set

Sr	Product No.	Name	Specification	QTY	Unit	Outview
----	-------------	------	---------------	-----	------	---------

1	09S1001	Pointed Reduction Forceps		1	pcs	
---	---------	---------------------------	--	---	-----	--



2	09S1002	Bent Reduction Forceps		1	pcs	
---	---------	------------------------	--	---	-----	--



3	09S1003	Centring Bone-holding Forceps		2	pcs	
---	---------	-------------------------------	--	---	-----	--



4	09S1004	Guide Wire	$\phi 1.2*150$	3	pcs
5	09S1005	Wire Sleeve	$\phi 1.2/ \phi 2.8$	1	pcs
6	09S1006	Threaded Drill Guide	$\phi 2.8$	3	pcs
7	09S1007	Drill Bit	$\phi 2.8*150$	2	pcs



8 09S1008 Limitation  $\phi$  2.8/SW2.5 2 pcs



9 09S1009 Wrench for Limitation SW2.5 1 pcs



10 09S1010 Plate Bender left 1 pcs



11 09S1011 Plate Bender right 1 pcs



12 09S1012 Screwdriver SW2.5 1 pcs



13 09S1013 Screwdriver SW2.5 1 pcs



14 09S1014 Sleeve 1 pcs



15 09S1015 Toque Limiting Screwdriver 1.5Nm/SW2.5 1 pcs



16 09S1016 Sleeve 1 pcs



17 09S1017 Alloy Drill  $\phi 4.0$  1 pcs



18 09S1018 Screwdriver for Removing the Screw SW2.5 1 pcs



19 09S1019 Drill Guide  $\phi 2.5$  1 pcs



20 09S1020 Drill Bit  $\phi 2.5$  1 pcs



21 09S1021 Drill Bit  $\phi 3.5$  1 pcs



22 09S1022 Drill Guide  $\phi 2.5 / \phi 3.5$  1 pcs



23 09S1023 Guide  $\phi 2.0$  1 pcs



24 09S1024 Guide  $\phi$  2.5 1 pcs



25 09S1025 Depth Gauge 1 pcs



26 09S1026 Tap HA3.5 1 pcs



27 09S1027 Tap HB4.0 1 pcs



28	09S1028	T-type Quick-coupling Handle		1	pcs
29	09S1029	Countersink Drill	Φ 6.0	1	pcs
30	09S1030	Extractor for Broken Screw	Φ 4.0	1	pcs
31	09S1031	Anti-slip extractor	SW2.5	1	pcs





32	09S1032	Anti-slip extractor		1	pcs
33	09S1033	Periosteal Detacher	round	1	pcs
34	09S1034	Periosteal Detacher	Flat	1	pcs
35	09S1035	Instrument Box		1	set
		In total		42	pcs



## LCP Small Fragment Instrument Set

