

# Proximal Femoral Nail Anti Rotation Instrument Set PFNA (PFNA II)

Sr	Product No.	Name	Specification	QTY	Unit	Outview
----	-------------	------	---------------	-----	------	---------

1	08S2001	Radiographic Ruler		1	pcs	 
2	08S2002	Protection Sleeve	$\phi 17 \times 120$	1	pcs	



3 08S2003 Pin Sleeve  $\phi 3.2$  1 pcs



4 08S2004 Proximal Drill  $\phi 16.5 / \phi 3.2$  1 pcs



5 08S2005 Quick Coupling Handle 1 pcs




6 08S2006 Quick Coupling 1 pcs



7 08S2007 Awl 1 pcs



8 08S2008 Reduction Rod 1 pcs




9 08S2009 Pin Holder 1 pcs



10 08S2010 Ball Tip Guide Pin  $\phi 2.5 / \phi 4.0 / 800$  1 pcs



11 08S2011 Reamer Tip  $\phi 8.5$  1 pcs



12 08S2012 Reamer Tip  $\phi 9$  1 pcs



13 08S2013 Reamer Tip  $\phi 9.5$  1 pcs



14 08S2014 Reamer Tip  $\phi 10$  1 pcs



15	08S2015	Reamer Tip	$\phi 10.5$	1	pcs
16	08S2016	Reamer Tip	$\phi 11$	1	pcs
17	08S2017	Reamer Tip	$\phi 11.5$	1	pcs



18 08S2018 Reamer Tip  $\phi$  12 1 pcs



19 08S2019 Reamer Tip  $\phi$  12.5 1 pcs



20 08S2020 Reamer Tip  $\phi$  13 1 pcs



21 08S2021 Reamer Tip Box  $\phi 8.5 - \phi 13$  1 pcs



22 08S2022 Reamer Shaft 2 pcs



23 08S2023 Insertion Handle 1 pcs





24	08S2024	Connection Screw	M12×1.5/SW7	2	pcs
25	08S2025	Wrench	SW7.0	1	pcs
26	08S2026	Connecting Rod	M12×1.5	1	pcs
27	08S2027	Combination Wrench	SW11	1	pcs



28	08S2028	Hammer		1	pcs
29	08S2029	Proximal Guide	180、200 Static	1	pcs
30	08S2030	Connection Screw	M8×1	3	pcs
31	08S2031	Sleeve for Blade	Φ 16/ Φ 11.2	1	pcs



32 08S2032 Compression Nut 1 pcs



33 08S2033 Sleeve  $\phi 11.2 / \phi 3.2$  1 pcs



34 08S2034 Pin  $\phi 3.2$  1 pcs



35	08S2035	Thead Pin	$\phi 3.2 \times 400$	5	pcs
36	08S2036	Depth Gauge	$\phi 3.2$	1	pcs
37	08S2037	Proximal Drill	$\phi 3.2 / \phi 10.5 / \phi 11.2$	1	pcs



38	08S2038	Proximal Drill	$\phi 10.5 / \phi 3.2$	1	pcs
39	08S2039	Limitation		1	pcs
40	08S2040	Screwdriver for Blade	SW4.5 / $\phi 3.2$	1	pcs



41 08S2041 Connecting Rod 1 pcs



42 08S2042 Distal Guide 1 pcs



43 08S2043 Sleeve  $\phi 11/\phi 8.2$  2 pcs



44	08S2044	Sleeve	$\phi 8.2 / \phi 4.0$	2	pcs
45	08S2045	Pin	$\phi 4.0$	1	pcs
46	08S2046	Guide Pin	$\phi 4.0 \times 300$	2	pcs
47	08S2047	Drill Bit	$\phi 4.0 \times 300$	3	pcs



48 08S2048 Limitation  $\phi$  4.0/SW3.5 1 pcs



49 08S2049 Depth Gauge 1 pcs



50 08S2050 Screwdriver SW3.5 1 pcs





51 08S2051 Wrench for End Cap SW7.0 1 pcs



52 08S2052 Flexible Screwdriver 1 pcs



53 08S2053 Blade Removal SW4.5/  $\phi$  3.2 1 pcs



54 08S2054 Nail Removal M12×1.5 M12×1.5 1 pcs



55 08S2055 Clean Pin  $\phi$  3.0 1 pcs



56 08S2056 Distal Connecting Rod for Long Nail 1 pcs



57	08S2057	Distal Guide for Long Nail		1	pcs
58	08S2058	Connecting Screw	M8×1	2	pcs
59	08S2059	Distal Position Device for Long Nail	left	1	pcs
60	08S2060	Distal Position Device for Long Nail	right	1	pcs



61	08S2061	Sleeve	$\phi 10 / \phi 8.1$	1	pcs
62	08S2062	Pin	$\phi 8.1$	1	pcs
63	08S2063	Sleeve	$\phi 8.1 / \phi 5.2$	1	pcs
64	08S2064	Drill Bit	$\phi 5.2$	1	pcs



65	08S2065	Drill Bit	$\phi 5.2$	1	pcs
66	08S2066	Position Rod	$\phi 8.1 / \phi 5.2$	1	pcs
67	08S2067	Position Block		1	pcs
68	08S2068	Instrument Box		1	set
		In total		82	pcs



## Proximal Femoral Nail Anti Rotation Instrument Set PFNA (PFNA II)

