




Universal Tibia Insterlocking Nail Instrument Set

Sr	Product No.	Name	Specification	QTY	Unit	Outview
1	08S3001	Radiographic Ruler		1	pcs	
2	08S3002	Thread Pin	$\phi 3.2 \times 300$	3	pcs	
3	08S3003	Protecting Sleeve	$\phi 17 \times 120$	1	pcs	

4 08S3004 Pin Sleeve $\phi 3.2$ 1 pcs



5 08S3005 Proximal Drill $\phi 12 / \phi 3.2$ 1 pcs



6 08S3006 Quick Coupling Handle 1 pcs



7 08S3007 Quick Coupling 1 pcs



8 08S3008 Awl 1 pcs



9 08S3009 Reduction Rod 1 pcs



10 08S3010 Pin Holder 1 pcs



11 08S3011 Ball Tip Pin $\phi 2.5 / \phi 4.0 / 800$ 1 pcs



12 08S3012 Reamer Tip $\phi 7.5$ 1 pcs



13 08S3013 Reamer Tip $\phi 8.0$ 1 pcs



14	08S3014	Reamer Tip	$\phi 8.5$	1	pcs
15	08S3015	Reamer Tip	$\phi 9$	1	pcs
16	08S3016	Reamer Tip	$\phi 9.5$	1	pcs
17	08S3017	Reamer Tip	$\phi 10$	1	pcs



18 08S3018 Reamer Tip ϕ 10.5 1 pcs



19 08S3019 Reamer Tip ϕ 11 1 pcs



20 08S3020 Reamer Tip ϕ 11.5 1 pcs



21 08S3021 Reamer Tip ϕ 12 1 pcs



22 08S3022 Reamer Tip Box $\phi 7.5 - \phi 12$ 1 pcs



23 08S3023 Flexible Reamer $\phi 7.2$ 2 pcs



24 08S3024 Insertion Handle 1 pcs



25 08S3025 Connection Screw M8×1/M6/SW6.5 2 pcs



26 08S3026 Wrench SW6.5 1 pcs



27 08S3027 Impact Block M8×1 1 pcs



28 08S3028 Nail Connector M8×1 2 pcs



29 08S3029 Wrench SW11 1 pcs



30 08S3030 Wrench $\phi 4.0$ 1 pcs



31 08S3031 Hammer Guide M8 \times 1/M10 \times 1 1 pcs



32 08S3032 Nut for Hammer Guide M10 \times 1 1 pcs



33 08S3033 Hammer 1 pcs



34 08S3034 Guide 1 pcs



35 08S3035 Connection Screw M8×1 1 pcs



36 08S3036 Position Sleeve $\phi 5.2$ 1 pcs



37 08S3037 Drill Bit $\phi 5.2$ 1 pcs



38 08S3038 Drill Bit $\phi 5.2$ 1 pcs



39 08S3039 Position Rod 1 pcs



40 08S3040 Position Block 1 pcs



41 08S3041 Distal Aiming Arm 1 pcs



42 08S3042 Connection Screw M11×1.5 1 pcs



43 08S3043 Sleeve 2 pcs



44 08S3044 Trocar 1 pcs



45 08S3045 Drill Sleeve $\phi 4.0$ 2 pcs



46 08S3046 Guide Pin $\phi 4.0 \times 300$ 2 pcs



47 08S3047 Distal Drill $\phi 4.0 \times 300$ 3 pcs



48 08S3048 Limitation $\phi 4.0 / SW3.5$ 1 pcs



49 08S3049 Temporary Position Rod ϕ 4.0 1 pcs



50 08S3050 Depth Gauge 1 pcs



51 08S3051 Screwdriver SW3.5 1 pcs



52 08S3052 Screwdriver SW3.5 1 pcs



53 08S3053 Proximal Aiming Arm 1 pcs



54 08S3054 Proximal Aiming Arm 1 pcs



55 08S3055 Connection Screw M6×1 2 pcs



56	08S3056	Sleeve	$\phi 11 / \phi 8.6$	1	pcs
57	08S3057	Sleeve	$\phi 8.6 / \phi 3.2$	1	pcs
58	08S3058	Trocar	$\phi 3.2$	1	pcs
59	08S3059	Drill Bit	$\phi 3.2$	2	pcs



60 08S3060 Limitation $\phi 3.2/SW3.5$ 1 pcs



61 08S3061 Compression Nut $\phi 4/M6/SW6.5$ 1 pcs



62 08S3062 Universal Wrench SW6.5 1 pcs



63 08S3063 Wrench for End Cap SW3.5 silicone handle 1 pcs



64	08S3064	End Cap Holder	SW3.5	1	pcs
65	08S3065	Instrument Box		1	set
		In total		77	pcs



Universal Tibia Insterlocking Nail Instrument Set

